# SAFETY DATA SHEET Xcem<sup>TM</sup> Septohex Instrument Cleanser & Disinfectant

<b>Document ID</b>	<b>Issue Date:</b>	
<u>Version</u>	Supersedes Date	

# > <u>IDENTIFICATION of SUBSTANCE, PREPARATION and COMPANY UNDERTAKING:</u>

1.1 Trade Name: Xcem<sup>TM</sup> Septohex Instrument Cleanser & Disinfectant

**1.2 Intended use of product:** For Instrument Cleaning & Disinfection.

**1.3. Xcem<sup>™</sup> Product ID:** This product can be purchased individually.

**1.4 Contact Address:** Xcem Ltd.

15 Picardy Road Belvedere Kent DA17 5QH

1.5 E-mail address and emergency contact number:

admin@xcem.co.uk

+44 (0) 203 490 7078

+44 (0) 203 490 2909

# **2. HAZARDS IDENTIFICATION:**

### Risk phrases:

Xcem<sup>™</sup> This product is not classified as hazardous according to EU Directive

# > 3. COMPOSITON/INFORMATION ON INGREDIENTS.

S.No	Ingredients	CAS.No.
1.0	Benzalkonium Chloride	68424-85-1
2.0	Didecyldimethylammonium chloride	7173-51-5
3.0	Chlorhexidine gluconate	55-56-1

# > 4. FIRST AID MEASURES:

- 4.1 Contact with eyes: Wash with large amount of water, if symptoms remains, get medical advice.
- 4.2 Contact with skin: Wash affected area with soap and water. If signs/symptoms persist and develop, get medical advice.

- 4.3 Inhalation: Remove person to fresh air. If signs/symptoms develop, get medical advice.
- 4.4 Ingestion: Do not induce vomiting unless instructed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention immediately.

# > 5. FIRE-FIGHTING MEASURES

- 5.1 Use fire extinguisher with class B extinguishing agent i.e. Dry chemical, Carbon dioxide.
- 5.2 Not flammable, ordinary combustible material.
- 5.3 Water may be used to blanket the fire. Use full Personnel Protective Equipment(PPE) and a self contained breathing apparatus.

# ➤ <u>6. ACCIDENTAL RELEASE MEASURES</u>

#### Personal Protection

- 1.0 Avoid the contact with the eyes
- 2.0 Don't let the material to come in contact with your skin
- 3.0 Don't touch the material
- 4.0 Don't ingest

# Method of Cleaning

- 1.0 Ventilate the area with fresh air
- 2.0 Clean any spillage
- 3.0 Cover with the absorbent.
- 4.0 Collect the maximum material as you can
- 5.0 Clean up residue
- 6.0 Dispose the collected material

#### > 7. STORAGE:

- 1.0 Store at Ambient temperature
- 2.0 Store away from Sun light
- 3.0 Store away from heat

# > 8 PHYSICAL AND CHEMICAL PROPERTIES

Appearance and Odour: whitish transparent gel

■ pH : 6.5-7..5 ■ Boiling Point : NA

Melting point : Not Determined
 Flash point : Not applicable
 Auto flammability : Not determined
 Water solubility : Not applicable
 Vapour Density : Not applicable
 Viscosity : Not determined

# > 9 STABILITY AND REACTIVITY

Stable at room temperature

# > 10 TOXICOLOGICAL INFORMATION

- **1.0** After contact with eyes moderate eyes irritation, symptoms may include redness, swelling, pain, tearing, and blurred
- 2.0 After contact with skin contact no significant irritation take place
- **3.0** After inhalation the signs/ symptoms may include cough, sneezing, nasal discharge, headache, hoarseness, nose and throat.
- 4.0 Sign of Ingestion include abdominal pain, stomach upset, nausea, Vomiting and diarrhea

### > 11 ECOLOGICAL INFORMATION

Environmental Data: Not determined.

Mobility in Soil and Water: Not determined. Persistence/Biodegradability: Not determined.

Ecotoxicity Data: Not determined.

# > 12 DISPOSAL CONSIDERATIONS

Dispose of wastes in sanitary landfill. Or incinerate waste product in an industrial or commercial incinerator.

The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Data sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications.