SAFETY DATA SHEET Xcem Ltd RESTOMER I

> 1. IDENTIFICATION of SUBSTANCE, PREPARATION and COMPANY UNDERTAKING:

1.1 Trade Name: Xcem Ltd's RESTOMER I

1.2 Intended use of product: Glass ionomer cement for Professional Use.

1.3. Xcem Ltd Product ID: 11352

1.4 Contact Address: Xcem Ltd.

15 Picardy Road Belvedere Kent DA17 5QH

1.5 E-mail address and emergency contact number:

admin@xcem.co.uk +44 (0) 203 490 7078 +44 (0) 203 490 2909

> 2. HAZARDS IDENTIFICATION:

Risk phrases: Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.

> <u>3. COMPOSITON/INFORMATION ON INGREDIENTS.</u>

Powder

S.No	Ingredients	CAS.No.
1.0	Glass Ionomer powder (Alumino Silicate Flouride)	
2.0	Poly Acrylic Acid (20%-30 %)	
3.0	Iron Oxide Pigment (Trace)	

> 4. FIRST AID MEASURES:

- 4.1 Contact with eyes: Wash with large amount of water, if symptoms remains, get medical advice.
- 4.2 Contact with skin: Wash affected area with soap and water. If signs/symptoms persist and develop, get medical advice.
- 4.3 Inhalation: Remove person to fresh air. If signs/symptoms develop, get medical advice.
- 4.4 Ingestion: Do not induce vomiting unless instructed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention immediately.

> 5. FIRE-FIGHTING MEASURES

- 5.1 Use fire extinguisher with class B extinguishing agent i.e. Dry chemical, Carbon dioxide.
- 5.2 Not flammable, ordinary combustible material.
- 5.3 Water may be used to blanket the fire. Use full Personnel Protective Equipment(PPE) and a self contained breathing apparatus.

> <u>6. ACCIDENTAL RELEASE MEASURES</u>

Personal Protection

- 1.0 Avoid the contact with the eyes
- 2.0 Don't let the material to come in contact with your skin
- 3.0 Don't touch the material
- 4.0 Don't ingest

Method of Cleaning

- 1.0 Ventilate the area with fresh air
- 2.0 Clean any spillage
- 3.0 Cover with the absorbent.
- 4.0 Collect the maximum material as you can
- 5.0 Clean up residue
- 6.0 Dispose the collected material

> 7. STORAGE:

- 1.0 Store at Ambient temperature
- 2.0 Store away from Sun light
- 3.0 Store away from heat

> 8 EXPOSURE CONTROL/PERSONAL PROTECTION

Occupational Exposure limits N/A Personal Protection PPE

Respiratory protection N/A

Skin Protection Overalls, rubber gloves
Eye protection Goggles/eye shields

Industrial Hygiene Wash hands thoroughly after contact with

this product

> 9 PHYSICAL AND CHEMICAL PROPERTIES

Appearance and Odour In different vita shades

Not Determined pН **Boiling Point** Not Determined Melting point Not Determined Flash point Not Determined Auto flammability Not Determined Water solubility Not Determined Vapour Density Not Determined Viscosity Not Determined

> 10 STABILITY AND REACTIVITY

Thermal decomposition / Conditions to avoid

Unstable at 98°C, gives acrylate vapours.

Product is stable in normal conditions.

Hazardous reactions / Substances to avoid

Avoid excessive heat and oxidising agents.

> 11 TOXICOLOGICAL INFORMATION

- 1.0 After contact with eyes moderate eyes irritation, symptoms may include redness, swelling, pain, tearing, and blurred
- 2.0 After contact with skin contact no significant irritation take place
- 3.0 After inhalation the signs/ symptoms may include cough, sneezing, nasal discharge, headache, hoarseness, nose and throat.
- 4.0 Sign of Ingestion include abdominal pain, stomach upset, nausea, Vomiting and diarrhea

> <u>12 ECOLOGICAL INFORMATION</u>

Environmental Data: Not determined.

Mobility in Soil and Water: Not determined. Persistence/Biodegradability: Not determined.

Ecotoxicity Data: Not determined.

> 13 DISPOSAL CONSIDERATIONS

Dispose of wastes in sanitary landfill. Or incinerate waste product in an industrial or commercial incinerator.

> 14 TRANSPORT INFORMATION

Powder:

UN-Number · ADR, ADN, IMDG,	not regulated	
IATA	C	
UN proper shipping name · ADR,	not regulated	
ADN, IMDG, IATA		
Transport hazard class(es) · ADR,	not regulated	
ADN, IMDG, IATA · Class		
Packing group - · ADR, IMDG,	not regulated	
IATA	not regulated	
Environmental hazards: · Marine	not regulated	
pollutant		
Special precautions for user	not regulated	
Transport in bulk according to		
Annex II of Marpol and the IBC	not regulated	
Code	C	
UN "Model Regulation"	not regulated	
O	\mathcal{C}	

> 15 REGULATORY INFORMATION

The information in this MSDS does not constitute an assessment of workplace risks.

> 16 OTHER INFORMALTION

The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Data sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications