SAFETY DATA SHEET Xcem Ltd TRYCEM NE

<u>1. IDENTIFICATION of SUBSTANCE, PREPARATION and COMPANY</u> <u>UNDERTAKING:</u>

1.1 Trade Name: Xcem Ltd's TRYCEM NE

1.2 Intended use of product	Zinc Oxide based temporary cement
	for Professional Use.

- **1.3. Xcem Ltd Product ID:** 4568
- **1.4 Contact Address:**Xcem Ltd.15 Picardy RoadBelvedere KentDA17 5QH

1.5 E-mail address and emergency contact number:

<u>admin@xcem.co.uk</u> +44 (0) 203 490 7078 +44 (0) 203 490 2909

> <u>2. HAZARDS IDENTIFICATION:</u>

Risk phrases: Unlikely to cause harmful effects under normal conditions of handling and use. Practice good hygiene. Wear gloves, overalls.

> <u>3. COMPOSITON/INFORMATION ON INGREDIENTS.</u>

Base Paste

S.No	Ingredients	CAS.No.
1.0	Zinc Oxide	1314-13-2

Catalyst Paste

S.No	Ingredients	CAS No
1.0	Organic Acids	

> <u>4. FIRST AID MEASURES:</u>

- 4.1 Contact with eyes: Wash with large amount of water, if symptoms remains, get medical advice.
- 4.2 Contact with skin: Wash affected area with soap and water. If signs/symptoms persist and develop, get medical advice.
- 4.3 Inhalation: Remove person to fresh air. If signs/symptoms develop, get medical advice.
- 4.4 Ingestion: Do not induce vomiting unless instructed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention immediately.

> <u>5. FIRE-FIGHTING MEASURES</u>

5.1 Use fire extinguisher with class B extinguishing agent i.e. Dry chemical, Carbon dioxide.

5.2 Not flammable, ordinary combustible material.

5.3 Water may be used to blanket the fire. Use full Personnel Protective Equipment(PPE) and a self contained breathing apparatus.

> <u>6. ACCIDENTAL RELEASE MEASURES</u>

Personal Protection

1.0 Avoid the contact with the eyes2.0 Don't let the material to come in contact with your skin3.0 Don't touch the material4.0 Don't ingest

Method of Cleaning

- 1.0 Ventilate the area with fresh air
- 2.0 Clean any spillage
- 3.0 Cover with the absorbent.
- 4.0 Collect the maximum material as you can
- 5.0 Clean up residue
- 6.0 Dispose the collected material

▶ <u>7. STORAGE:</u>

- 1.0 Store at Ambient temperature
- 2.0 Store away from Sun light
- 3.0 Store away from heat

> <u>8 EXPOSURE CONTROL/PERSONAL PROTECTION</u>

Occupational Exposure limits	N/A
Personal Protection	PPE
Respiratory protection	N/A
Skin Protection	Overalls, rubber gloves
Eye protection	Goggles/eye shields
Industrial Hygiene	Wash hands thoroughly after contact with this product

> 9 PHYSICAL AND CHEMICAL PROPERTIES

Appearance and Odour	Catalyst Light Brown Coloured paste with characteristic smell. Base White paste with no specific smell
pН	Not Determined
Boiling Point	Not Determined
Melting point	Not Determined
Flash point	Not Determined
Auto flammability	Not Determined
Water solubility	Not Determined

Not Determined

Not Determined

> 10 STABILITY AND REACTIVITY

Vapour Density

Viscosity

Reactivity: No specific reactivity is known. Chemical stability · Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications. Possibility of hazardous reactions: No dangerous reactions known Conditions to avoid: No further relevant information available. Incompatible materials: No further relevant information available. Hazardous decomposition products: No dangerous decomposition products known

> <u>11 TOXICOLOGICAL INFORMATION</u>

- 1.0 After contact with eyes moderate eyes irritation, symptoms may include redness, swelling, pain, tearing, and blurred
- 2.0 After contact with skin contact no significant irritation take place
- 3.0 After inhalation the signs/ symptoms may include cough, sneezing, nasal discharge, headache, hoarseness, nose and throat.
- 4.0 Sign of Ingestion include abdominal pain, stomach upset, nausea,

Vomiting and diarrhea

> <u>12 ECOLOGICAL INFORMATION</u>

Environmental Data: Not determined. Mobility in Soil and Water: Not determined. Persistence/Biodegradability: Not determined. Ecotoxicity Data: Not determined.

> <u>13 DISPOSAL CONSIDERATIONS</u>

Dispose of wastes in sanitary landfill. Or incinerate waste product in an industrial or commercial incinerator.

> <u>14 TRANSPORT INFORMATION</u>

UN-Number · ADR, ADN, IMDG, IATA UN proper shipping name · ADR, ADN, IMDG, IATA Transport hazard class(es) · ADR, ADN, IMDG, IATA · Class Packing group - · ADR, IMDG, IATA Environmental hazards: · Marine pollutant Special precautions for user Transport in bulk according to Annex II of Marpol and the IBC Code UN "Model Regulation"

> <u>15 REGULATORY INFORMATION</u>

The information in this MSDS does not constitute an assessment of workplace risks.

> <u>16 OTHER INFORMALTION</u>

The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Data sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications.