SAFETY DATA SHEET Xcem Ltd.'s Imprasil HB

<u>1. IDENTIFICATION of SUBSTANCE, PREPARATION and COMPANY</u> <u>UNDERTAKING:</u>

1.1 Trade Name: Xcem Ltd's Imprasil HB

1.2 Intended use of product: Dental impression material Professional Use.

1.3. Xcem Ltd Product ID:

1.4 Contact Address:

Xcem Ltd. 15 Picardy Road Belvedere Kent DA17 5QH

1.5 E-mail address and emergency contact number:

<u>admin@xcem.co.uk</u> +44 (0) 203 490 7078 +44 (0) 203 490 2909

> <u>2. HAZARDS IDENTIFICATION:</u>

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP] Not classified

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Extra phrases: The product is seen as a medical device and therefore not subject to labelling (EU- regulation 1907/2006, article 2, paragraph 6c).

2.3 Other Hazards

Other hazards not contributing to the : None under normal conditions. classification

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

> <u>3. COMPOSITON/INFORMATION ON INGREDIENTS.</u>

Ingredients	CAS.No.	Percentage
Hydroxyl terminated silicone		20%-50%
High Purity Silica		40%-50%

> <u>4. FIRST AID MEASURES:</u>

- 4.1 Contact with eyes: Wash with large amount of water, if symptoms remains, get medical advice.
- 4.2 Contact with skin: Wash affected area with soap and water. If signs/symptoms persist and develop, get medical advice.
- 4.3 Inhalation: Remove person to fresh air. If signs/symptoms develop, get medical advice.
- 4.4 Ingestion: Do not induce vomiting unless instructed to do so by medical personnel. Never give anything by mouth to an unconscious person. Get medical attention immediately.

> <u>5. FIRE-FIGHTING MEASURES</u>

5.1 Use fire extinguisher with class B extinguishing agent i.e. Dry chemical, Carbon dioxide.

5.2 Not flammable, ordinary combustible material.

5.3 Water may be used to blanket the fire. Use full Personnel Protective Equipment(PPE) and a self contained breathing apparatus.

> <u>6. ACCIDENTAL RELEASE MEASURES</u>

Personal Protection

- 1.0 Avoid the contact with the eyes
- 2.0 Don't let the material to come in contact with your skin
- 3.0 Don't touch the material
- 4.0 Don't ingest

Method of Cleaning

- 1.0 Ventilate the area with fresh air
- 2.0 Clean any spillage
- 3.0 Cover with the absorbent.
- 4.0 Collect the maximum material as you can
- 5.0 Clean up residue
- 6.0 Dispose the collected material

> <u>7. STORAGE:</u>

1.0 Store at Ambient temperature

- 2.0 Store away from Sun light
- 3.0 Store away from heat

> <u>8 EXPOSURE CONTROL/PERSONAL PROTECTION</u>

Occupational Exposure limits N/A Personal Protection PPE

Respiratory protection Skin Protection Eye protection Industrial Hygiene N/A Overalls, rubber gloves Goggles/eye shields Wash hands thoroughly after contact with this product

> <u>9 PHYSICAL AND CHEMICAL PROPERTIES</u>

Appearance and Odour pH Boiling Point Melting point Flash point Auto flammability Water solubility Vapour Density Viscosity Solid Putty, No characteristic smell Not Determined Not Determined

> <u>10 STABILITY AND REACTIVITY</u>

Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

Chemical stability

Stable under normal conditions.

Possibility of hazardous reactions No polymerization.

Conditions to avoid No flames, no sparks. Eliminate all sources of ignition.

Incompatible materials Oxidizing agent.

Hazardous decomposition products No decomposition if stored normally.

> <u>11 TOXICOLOGICAL INFORMATION</u>

- 1.0 After contact with eyes moderate eyes irritation, symptoms may include redness, swelling, pain, tearing, and blurred
- 2.0 After contact with skin contact no significant irritation take place
- 3.0 After inhalation the signs/ symptoms may include cough, sneezing, nasal discharge, headache, hoarseness, nose and throat.
- 4.0 Sign of Ingestion include abdominal pain, stomach upset, nausea, Vomiting and diarrhea

> <u>12 ECOLOGICAL INFORMATION</u>

Environmental Data: Not determined. Mobility in Soil and Water: Not determined. Persistence/Biodegradability: Not determined. Ecotoxicity Data: Not determined.

> <u>13 DISPOSAL CONSIDERATIONS</u>

Dispose of wastes in sanitary landfill. Or incinerate waste product in an industrial or commercial incinerator.

> <u>14 TRANSPORT INFORMATION</u>

In accordance with ADR / RID / IMDG / IATA / ADN UN number Not regulated for transport UN proper shipping name Transport hazard class (es) Packing group Environmental hazards Dangerous for the environment: No Other information: No supplementary information available Special precautions for user Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code Not applicable

> <u>15 REGULATORY INFORMATION</u>

The information in this MSDS does not constitute an assessment of workplace risks.

> <u>16 OTHER INFORMALTION</u>

The information on this Safety Data Sheet is based on our experience and is correct to the best of our knowledge at the date of publication, but we do not accept any liability for any loss, damage or injury resulting from its use (except as required by law). The information may not be valid for any use not referred to in this Data sheet or use of the product in combination with other materials. For these reasons, it is important that customers carry out their own test to satisfy themselves as to the suitability of the product for their own intended applications.